

# PREPCLEAN®

## Forever Wet

### Instrument Transport Humectant Pre-Cleaner

#### Non-Enzymatic, Non-Aerosol

#### Helps to Prevent the Adhesion of Bio-burden to Surgical Instruments and Scopes

Prepclean® Forever Wet is a neutral pH, non-aerosol humectant spray which promotes the long lasting retention of moisture on soiled instruments helping to prevent the adhesion of bio-burden to surgical instruments and scopes. This unique formulation gently coats instruments to maintain moisture making it an ideal pre-cleaner for soiled instruments during transport or when left for an extended period of time.

Perfect for use in Operating Theatres, Endoscopy Suites, Outpatient Surgery, Dental and other departments where instruments are transported to decontamination.

- Humectant spray promotes long lasting retention of moisture on soiled instruments helping to eliminate dried on bio-burden
- Keeps instruments moist for an extended period of time
- Perfect for instruments that set overnight or over a weekend
- No "skinning"; unlike other pre-cleaners does not leave a hardened skin on surfaces after an extended period of time
- Helps to prevent the adhesion of blood clots and organic debris in the cannula tunnel of instruments and scopes
- Begins working instantly on contact
- Rinses easily; washer/disinfector compatible
- No messy spills or splashing during transportation
- Neutral pH; safe for all instruments and scopes
- Non-aerosol; environmentally safe

#### Available in:

Cat No 345PCF32 6 x 946ml



#### DIRECTIONS FOR USE

1. Spray directly on the soiled instruments and scopes immediately after use.
2. Allow the instrument to sit in the Prepclean until ready for cleaning. As long as all surfaces have been completely covered it is not necessary to reapply as moisture is retained.
3. Rinse instruments before cleaning, unless they are being cleaned in a Washer Disinfector.

#### CAUTIONS:

Avoid prolonged contact with skin.  
Do not swallow. Keep away from children.