



FINAL REPORT

CLEANING EVALUATION VALIDATION

PROCEDURE NO. STP0083 REV 04
PROTOCOL DETAIL SHEET NO. 200900083 REV 01

LABORATORY NO. 462879

PREPARED FOR:

FRANK BASS
RUHOF HEALTHCARE
393 SAGAMORE AVE.
MINEOLA NY 11501

SUBMITTED BY:

NELSON LABORATORIES, INC.
6280 S. REDWOOD RD.
SALT LAKE CITY UT 84123-6600
801-290-7500



CLEANING EVALUATION VALIDATION

LABORATORY NUMBER:	462879
PROCEDURE NUMBER:	STP0083 REV 04
PROTOCOL DETAIL SHEET NUMBER:	200900083 REV 01
SAMPLE SOURCE:	Ruhof Healthcare
SAMPLE IDENTIFICATION:	(Endozime XP Batch Number 7186 Manufactured January 22, 2009 (1 bottle)) Sponsor Supplied, 1 box Proformance (TOSI) automated washer test coupons, Sponsor Supplied P.O. #Ruhof101
DEVIATIONS:	None
NUMBER OF TEST SAMPLES:	30 (Nelson Laboratories, Inc. provided 4 additional coupons)
PROTOCOL APPROVAL DATE:	11 Feb 2009
SAMPLE RECEIVED DATE:	17 Feb 2009
LAB PHASE START DATE:	24 Feb 2009
LAB PHASE COMPLETION DATE:	02 Mar 2009
REPORT ISSUE DATE:	03 Mar 2009

INTRODUCTION:

A cleaning validation was performed using the manufacturer's recommended cleaning procedures and specifications. One positive coupon and three test coupons were utilized to determine the reduction in bioburden. An extraction was performed on the positive coupon. Three sponsor provided TOSI coupons were also cleaned and visually inspected.

PROCEDURES:

Culture Preparation: The defibrinated blood soil (DBLSO) was inoculated with *Geobacillus stearothermophilus*, ATCC #7953 from a stock spore suspension maintained at 2-8°C to yield a minimum population of 10^4 colony forming units (CFU)/mL. A standard plate count was performed on the inoculated test soil to determine the initial titer of the test organism.


Sample Contamination: The coupons were immersed in inoculated soil and allowed to remain in contact with the soil for a minimum of 15 minutes. The soiled coupons were placed into a clean pan and the pan was covered with a towel dampened with sterile purified water (PURW). The coupons were allowed to set for a minimum of 30 minutes to simulate the wait time between contamination and cleaning.

Ruhof Healthcare
Lab Number 462879

Cleaning Evaluation Validation

Water Hardness Test: The water hardness was measured the day of automated procedure testing. A sample of cold water from the water inlet and a sample of hot water from the water inlet were collected in separate sterile containers. The water hardness was measured using a Model 5B Hardness Test Kit.

Cleaning Procedure: The coupons were transferred into mesh baskets and loaded onto the 4-Level manifold rack accessory contained inside the STERIS® 444 washer for processing. Sample placement as well as the dunnage placement used to simulate a maximum load is shown in the following figure:



Dunnage = ~20 scissors

The cycle was selected and the following set of parameters was properly programmed:

PHASE	RECIRCULATION TIME (MINUTES)	WATER TEMPERATURE	DETERGENT TYPE AND CONCENTRATION
Pre-Wash 1	03:00	Cold tap water	N/A
Enzyme Wash	03:00	Hot tap water	Detergent provided by the sponsor (1/4 oz/gal)
Wash 1	03:00	65.5°C (Set Point)	Detergent provided by the sponsor (1/4 oz/gal)
Rinse 1	03:00	Hot tap water	N/A
Thermal Disinfection	05:00	90°C	N/A

Upon completion of the draining phase after Rinse 1, the cycle was stopped and the coupons were removed from the washer.

The coupons were visually examined by the naked eye under normal lighting conditions to determine if all adherent visible soil (e.g. blood, protein substances and other debris) had been removed from the surfaces, lumens, crevices, and serrations. Refer to Figure 1.

Bioburden Testing: The positive and cleaned coupons were tested for bioburden by immersing in peptone Tween® (PEPT) and shaking manually 100 times to extract the organisms present. Aliquots of the extract fluid were diluted where appropriate, and plated onto soybean casein digest agar (SCDA) or filtered through a 0.45 µm membrane and the membrane was placed onto SCDA. Plates were incubated at 55-60°C for 24 ± 3 hours and colonies were enumerated.

Calculations for the Bioburden: The percent reduction in bacterial counts was calculated using the following formula:

$$\% \text{ Reduction} = 100 - \left(\frac{\text{final population}}{\text{initial population}} \times 100 \right)$$

Log reductions were calculated using this formula:

$$\text{Log Reduction} = \text{Log Initial Population} - \text{Log Final Population}$$

Where: Final population = Number of organisms recovered from the cleaned coupons
Initial population = Positive coupon titer

Estimations were reported for counts which were outside of the statistically accurate range.

Ruhof Healthcare
Lab Number 462879

Cleaning Evaluation Validation

RESULTS:

The results of the cleaning evaluation are summarized in Table 1. Each coupon was individually inoculated, and some variability is expected between coupons. The soil titer and the positive coupon titer are listed at the bottom of the table. The TOSI coupons were compared to the TOSI chart, and all coupons were a 0 according to the TOSI chart listed in Figure 1.

Water Hardness Testing Results: Cold Water was 222.3 ppm
Hot Water was 188.1 ppm

STATEMENT OF UNCERTAINTY:

If applicable, the statement of uncertainty is available to sponsors upon request.



Technical Reviewer



Nick Workman
Study Director



Study Completion Date

emm

0-5 T.O.S.I. chart scale




Figure 1. T.O.S.I. Chart

TABLE 1. Cleaning Results

DEVICE NUMBER	COUNTS PER DEVICE (CFU)	PERCENT REDUCTION (%)	LOG ₁₀ REDUCTION
1	1.1×10^0	99.99979	5.7
2	$<1.1 \times 10^0$	>99.99979	>5.7
3	$<1.1 \times 10^0$	>99.99979	>5.7

Soil Titer: 1.0×10^4 CFU/mL

Positive Device Titer: 5.2×10^5 CFU/device



Ruhof Healthcare
Lab Number 462879

Cleaning Evaluation Validation

All reports and letters issued by Nelson Laboratories, Inc. are for the exclusive use of the sponsor to whom they are addressed. These results relate only to the samples tested. Reports may not be reproduced except in their entirety. No quotations from reports or use of the corporate name is permitted except as expressly authorized by Nelson Laboratories, Inc. in writing. The significance of any data is subject to the adequacy and representative character of the samples tendered for testing. Nelson Laboratories, Inc. warrants that all tests are performed in accordance with established laboratory procedures and standards. Nelson Laboratories, Inc. makes no other warranties of any kind, express or implied. Nelson Laboratories, Inc. expressly states that it makes no representation or warranty regarding the adequacy of the samples tendered for testing for any specific use of application, that determination being the sole responsibility of the sponsor. Nelson Laboratories' liability for any loss or damage resulting from its actions or failure to act shall not exceed the cost of tests performed, and it shall not be liable for any incidental or consequential damages.