

# BIO-CLEAN

## Phase One

### NON-ENZYMATIC NEUTRAL DETERGENT ENDOSCOPY BEDSIDE CARE KIT

The only Bedside Care Kit that removes gross contaminants, including synthetic lipids while pre-cleaning scope surface and internal channels enabling high level disinfectant to kill and remove biofilm



### Scope Cleaning Begins at the Bedside!

The kit is designed for convenient bedside use to remove all contaminants from the outside sheath and inside channel of scopes. The kit contains an Bio-Clean Sponge to remove soil from the insertion tube and Bio-Clean Non-Enzymatic neutral detergent solution for suctioning through the channels. The kit provides optimal cleaning, preventing dried on soil, before the scope is transported for reprocessing. The Bio-Clean Phase One with Bio-Clean Technology is designed to remove synthetic lipids while targeting insoluble polysaccharides that encase Biofilm allowing for the complete elimination of all Bioburden and Biofilm by high level disinfectants.

### Features and Benefits:

- Begins cleaning on contact, preventing soil from drying on scope surfaces.
- Removes gross contaminants, including synthetic lipids.
- Uniquely fitted for all scopes: flexible/rigid.
- Safe on all scopes.
- Neutral pH, non-abrasive, 100% Biodegradable.
- Empty containers can be rinsed and recycled.
- Will not harm any metals, plastic, rubber, corrugated tubing, glass or mirrors.
- Each kit is individually packaged for one time use, eliminating the risk of cross-contamination.

### Contents:

Bio-Clean Sponge - contoured sponge for cleaning outside sheath of scope.

Bio-Clean Non-Enzymatic Neutral Detergent - 250ml pre-diluted.

### DIRECTIONS FOR USE:

On the container label.

### Available in:

Cat No 345BPO250 250ml

Unit of issue - 24