

ENDOZIME SPONGE®

Contoured sponge for cleaning purposes. A great leap forward in the cleaning and maintaining of your Endoscopes.

These sponges are designed to contour the outside sheath of a scope and are soaked in Endozime to begin cleaning on contact.

Suitable for removing gross contaminants before enzymatic pre-soak.

Endozime sponges are designed to provide a quicker and more thorough cleaning of endoscopes and flexible fibreoptic or tubular medical instruments.



AVAILABLE IN

Cat No 345SPG
Unit of issue - 100



Here are some of the features and benefits of this sponge:

Easy to use - specifically contoured to fit endoscopes and other tubular instruments.

Convenient and safe – each Endozime Sponge is packed individually for one-time use.

Uniquely designed to thoroughly clean all surface areas of rigid or flexible endoscopes including cysto-urethrosopes, laparoscopes, proctoscopes, etc.

Hospital-proven enzymatic detergent action – safely and effectively removes organic and proteinaceous materials, significantly reducing micro-organisms.

Increases effectiveness of manual scope cleaning – because speed and simplicity of use allows you to remove gross contaminants while scopes are still wet.

Does not adversely affect delicate endoscopic parts, such as lenses, lighting fixtures or air/water nozzles.

100% biodegradable.

Eliminates biological odours on contact.

R **RUHOF** 
THE EXPERTS IN A BOTTLE

Peskett Solutions Ltd
Unit 4-5 Brampton Business Park
55 Brampton Road, Eastbourne, East Sussex, BN22 9AF
Tel: 01323 511038 Fax: 01323 503247
E Mail: info@ruhof.co.uk Web: www.ruhof.co.uk