

RINSE AID

DESIGNED FOR USE IN WASHER-DISINFECTORS, TUNNEL WASHERS AND CART WASHERS.

Rinse Aid has not been tested for being free-rinsing, but the purpose of Rinse Aid is exactly that to further lower the surface tension and thus provide a surface free of cleaning product residue.

Description

For surgical instruments, polyware and glassware.

Contains powerful wetting agents, which lower the surface tension of the rinse solution.

Rinse water drains quickly and completely.

Prevents water spotting.

Leaves surfaces clean and sparkling.

The answer to your hard water problems.

DIRECTIONS FOR USE

- Add 1 to 3ml per Litre of water in the final rinse of each automatic wash load (dependent on water quality).

Available in:

Cat No 345RN	4 x 4 Litres
Cat No 345RN10L	2 x 10 Litres



Peskett Solutions Ltd
Unit 4-5 Brampton Business Park
55 Brampton Road, Eastbourne, East Sussex, BN22 9AF
Tel: 01323 511038 Fax: 01323 503247
E Mail: info@ruhof.co.uk Web: www.ruhof.co.uk